Cleaning protocol

How to clean the Simple Up. Also follow your organisation's guidelines for cleaning and infection prevention.

✓ Approved agents

For cleaning the surface, Melius Pro approves products containing diluted solutions of:

- cocoalkyl dimethylbenzyl ammonium chloride;
- chloroxylenol.

These substances are found in, for example, Lysol, Dettol and Tuffie wipes.

Also suitable: mild soap, such as washing-up liquid.

✗ Do not use

The Simple Up cannot withstand:

- disinfection in a steriliser or autoclave;
- agents or devices that disinfect using heat or steam;
- immersion in water, of the device or of any parts.

When and what

When	What
Between patients	Wipe the grab handles and the standing surface, preferably with disinfectant wipes.
Monthly, or when visibly dirty	Rubber anti-slip mat: remove loose dirt with a soft brush and wipe with a damp cloth.
Regularly, and at least once a year	Clean the housing, preferably with disinfectant wipes. Remove dust from the bellows and wipe with a damp cloth.

For inspection and maintenance, see the user manual. If the anti-slip mat is worn or damaged, have it replaced (part number CL-015).

Issued by Melius Pro B.V., manufacturer of the Simple Up.

Replaces the cleaning protocol dated 29 March 2013. The approved agents are unchanged.

Melius Pro B.V.
Hamersveldseweg 135a, 3833 GN Leusden, the Netherlands
+31 33 434 50 60 · info@meliuspro.com · www.meliuspro.com